



Double antibacterial effects of silver and copper for long-lasting clean use

Anti-Bacterial Bottle Sponge



introduction

With the antibacterial effects of silver and copper, it uses for hygienic and reach the edge of the bottle to remove dirt.



1
Long-lasting hygienic use with double antibacterial effects of silver and copper

The double antibacterial effects of silver and copper control of bacteria on sponges.



2
Hygienic because the handle is also treated with an antibacterial coating

The silver antimicrobial agent in the handle makes it possible to use not only the sponge but also the handle clean.



3
Rectangular shape to reach the edge of the bottle

The shape of the bottle allows for thorough cleaning to the very tip of the bottle, and it can be installed and used in either direction.



4
Good foaming and drainage

A small amount of detergent produces plenty of lather, and the sponge is easy to drain, so it can be used clean.



5
Replacement sponge for long-lasting and economical


A replacement sponge is also available, so only the sponge can be replaced when it gets dirty, making it economical and long-lasting.



The position of the sponge can be adjusted and oriented to suit the application, and the length of the sponge can be adjusted for easy washing.



Firmly clamp the sponge with the sponge fixing claws and slide the clasp to the sponge side to fix it



To make it shorter, fix the sponge at the back of the handle; to make it longer, fix it at the end of the handle.



The sponge can be installed in either orientation.

